

INTEROSS[®] COLLAGEN

ANORGANIC BONE-COLLAGEN COMPOSITE



WHAT IS INTEROSS[®] COLLAGEN

InterOss[®] Collagen is an anorganic hydroxyapatite-collagen composite for use in periodontal, oral and maxillofacial surgery. It is a combination of bovine bone granules and collagen fibers for exceptional handling. Its interconnected porous structure and osteoconductive property allows for the ingrowth of the adjacent viable bone while permitting easy exchange of body fluids. InterOss[®] Collagen is available in various sizes in plug and block forms.



InterOss[®] Collagen Block-Type

Reference	IOC-50	IOC-100	IOC-250	IOC-350	IOC-500
Weight	50mg	100mg	250mg	350mg	500mg



SIGMAGRAFT

THE BONE GRAFT SUBSTITUTE SPECIALISTS

INDICATIONS FOR USE


InterOss® Collagen is indicated for the filling of extraction sockets to enhance the preservation of the alveolar ridge. The product is recommended for:

- ✔ Filling of extraction sockets to enhance preservation of the alveolar ridge
- ✔ Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) & Guided Bone Regeneration (GBR)

BENEFITS



Formable and spongy



Highly porous



Easily trimmable

- ✔ Outstanding porosity allows for rapid bone regeneration and integration to the host tissue
- ✔ Exhibit excellent mechanical property providing enhanced support to the newly forming bone
- ✔ High water absorption capability allows for rapid blood infiltration and nutrient exchange
- ✔ Spongy nature allows for ease of trimming and handling

Contact Us

Address: SigmaGraft, Inc. 575 Sally Place Fullerton, CA 92831, USA
Phone: (714) 525-0114, Toll Free: (888) 499-0114
E-Mail: info@sigmagraft.com